

## Microprocessor pH Conductivity



**Product Categories:** [Analytical Lab Instruments](#)

**Product Tags:** [Analytical Laboratory Equipment](#)

**Product Page:**

<https://www.labappara.com/product/microprocessor-ph-conductivity/>

### Product Description

Microprocessor pH Conductivity - Analytical Laboratory Equipment

Measurement of pH, conductivity-mV & temperature 16x2 Alphanumeric LCD

Display with pH combination electrode, conductivity cell and temperature probe.

PC Interface facility.

#### **Accessories On Request**

- a) pH combination electrode
- b) Gel filled epoxy body combination electrode
- c) Glass electrode
- d) Silver electrode
- e) Platinum electrode
- f) Reference electrode
- g) pH electrode stand
- h) Temperature probe
- i) Duracell battery or equivalent (9V DC)

j) ORP Electrode

k) Standard pH Buffer Solution