

Microprocessor Spectrophotometer



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Laboratory Equipment](#)

Product Page:

<https://www.labappara.com/product/microprocessor-spectrophotometer-2/>

Product Description

Microprocessor Spectrophotometer - Analytical Laboratory Equipment
340-960 nm with facility for automatic concentration, % transmission and absorbance measurements.

16x2 alphanumeric LCD display for wavelength resolution.

1 nm Cuvettes holder provision for 10 mm and 50 mm path length cuvettes with a set of 1 pair.

Visible range 340-1000 nm with facility for automatic concentration, % Transmission and Absorbance measurement.

Wavelength resolution 1 nm 3 and half digit LCD display.

320 -1000 nm & 16x2 line backlit LCD Display with PC Attachment Facility Storage for 1000 samples, RS 232C Computer Interface.

a) Glass Cuvettes 10 mm path length (set of 2)

- b) Glass Cuvette 50 mm (set of 2)
- c) Glass Cuvette 100 mm (set of 2)
- d) Cuvette holder 100 mm path length
- e) Cuvette holder 50 mm path length
- f) Spare lamp (Vis.) Halogen 6V 10W
- g) Polystyrene Cuvettes 10mm path length (set of 5)