

MICROPROCESSOR SPECTROPHOTOMETER



Product Categories: [Analytical Lab Instruments](#), [Spectrophotometer](#)

Product Tags: [Analytical Laboratory Equipment](#), [Spectrophotometer](#)

Product Page: <https://www.labappara.com/product/microprocessor-spectrophotometer/>

Product Description

Microprocessor Spectrophotometer - Analytical Laboratory Equipment

340-960 nm with facility for automatic concentration, % transmission and absorbency measurements.

16x2 alphanumeric LCD display for wavelength resolution 1 nm

Cuvettes holder provision for 10 mm and 50 mm path length cuvettes with a set of 1 pair.

We are counted amidst the well-known Manufacturers, Exporters & Suppliers of highly efficient Microprocessor Spectrophotometer.

It is a solid state micro controller based instrument that is specifically designed to provide quick and accurate analysis of any concentration.

By using the latest microprocessor technology and advanced engineering techniques, our Spectrophotometer works with in the wavelength ranging from 320 nm to 1100 nm.