

# Microprocessor Turbidity Nephelo Meter



**Product Categories:** [Analytical Lab Instruments](#)

**Product Tags:** [Analytical Laboratory Equipment](#)

**Product Page:**

<https://www.labappara.com/product/microprocessor-turbidity-nephelo-meter/>

## Product Description

Microprocessor Turbidity Nephelo Meter - Analytical Laboratory Equipment

Four digit LED display Range up to 1000 NTU/JTU.

Resolution 0.1 NTU/JTU.

Soft touch membrane keys operated.

a) Test tube set of 5 850

b) Spare lamp 6.8V, 0.3 Amp. / Source LED 450.Auto ranging facility.

c) Standard Solution Turbidity 400 NTU (100 ml) .

Display

4 Digit LED Display

Detector

Photocell/ Photodiode

## Range

0 to 200, 1000 NTU/JTU

## Resolution

1 NTU/JTU

## Auto Standard Calibration

5 Point Calibration

## Accuracy

1% FS 1 Digit

## Sample System

30 mm Clear glass test tube

## Power Supply

230V 10%, 50HZ AC