

MICROPROCESSOR VISCOMETER



Product Categories: [Analytical Lab Instruments](#)

Product Tags: [Analytical Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/microprocessor-viscometer-2/>

Product Description

MICROPROCESSOR VISCOMETER - Analytical Laboratory Equipment

Having advantage of high accuracy, stable digital display easy to operate and good anti-interference.

It is widely used in the measurement of viscosity of various kinds of fluids such as grease, paints, food, medicine, adhesive and cosmetic etc.

It can also be used for liquids.

LCD display measurement ranges 1,00,000 mPa.s with 4 spindles.