

Microscope Lamp (LED)



Product Categories: [Microscope](#), [Microscope Assecories](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/microscope-lamp-led-type/>

Product Description

Microscope Lamp (LED) - Scientific Equipment

LED emit light with a relatively narrow wavelength band, white light can be obtained by combining an LED and a fluorescent material or a combination of multiple LED . It features compact size, low power consumption, and long life, but the power may be poor at specific wavelengths depending on how the LED are combined.