

Microtip Box



Product Categories: [Laboratory Plastic Ware](#), [Vails](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/microtip-box/>

Product Description

Microtip Box - Scientific Laboratory Equipment

Micro Tip Box: These unique boxes are moulded in polypropylene and are therefore autoclavable. The box itself is clear providing an excellent view of the tips inside while a table like structure is provided inside to hold the tips. The tip holding spaces are uniformly placed to facilitate the use of multi channel pipettes. The cover is joined to the base through built in hinge and provides excellent protection to the tips.

Features:

Flawless finish

Precisely designed

Durability

Particulars

Packing

For 96 Micro Tips of 2-10 μ l

6 Pcs

For 96 Micro Tips of 2-200 μ l

6 Pcs

For 100 Micro Tips of 200-1000 μ l

6 Pcs