

Milli Ampere Meter M.I. (Single Range) 500mA



Product Categories: [Electronics](#), [Engineering Equipment](#), [Moving Iron AC/DC Voltmeters And Ampere meters](#), [Moving Iron AC/DC Voltmeters And Ampere meters](#)

Product Page:

<https://www.labappara.com/product/milli-ampere-meter-m-i-single-range-500ma/>

Product Description

Milli Ampere Meter M.I. (Single Range) 500mA

Specification:

High quality bench type instrument for the use schools and education

Durable ABS plastic case

No power supply required for measurement

With mirror scale and 4 mm safety jack. Class 2,0

Range 0-50mA / 500mA / 5A DC

Accuracy $\pm 2,0$ % of fullscale

Dimensions: 90x106x103mm