

Millman's Theorem



Product Categories: [Electronics](#), [Engineering Equipment](#), [Network Theorem's Verifications \(In DC Circuits\)](#), [Network Theorem's Verifications \(In DC Circuits\)](#)

Product Page: <https://www.labappara.com/product/millmans-theorem/>

Product Description

Millman's Theorem

Specifications :

Instrument comprises of three fixed output DC regulated power supplies of 12V, 15V & 18V. Two moving coil analog meter are mounted on the front panel to measure voltage & current, circuit diagram is printed & important connections are brought out on sockets.