

# Mirror Strip



**Product Categories:** [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/mirror-strip/>

## Product Description

Mirror Strip - Physics Equipment

Mirror Strip are specifically process as per required sizes from true plain Mirror, having a silver coating with the protective layer of copper and led with the additional layer of epoxy paint, used in scientific application for the reflection of light. Silver back with protective coating.

Size:-

(a) 75 x 50 m

(b) 100 x 50 mm

(c) 100 x 25 mm