

Model Of Terrestrial Telescope (Fitted On Base)



Product Categories: [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/model-terrestrial-telescope-fitted-base/>

Product Description

Model Of Terrestrial Telescope (Fitted On Base) - Physics Optical Equipment

An optical telescope which uses lenses is known as a refracting telescope or a Refractor one which uses a mirror is known as a reflecting telescope or a Reflector. Besides optical telescopes, astronomers also use telescopes that focus radio waves, X-rays, and other forms of electromagnetic radiation.