

Mosfet Amplifier



Product Categories: [Amplifier's](#), [Amplifier's](#), [Electronics](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/mosfet-amplifier/>

Product Description

Mosfet Amplifier

Features:

Instrument comprises of two dc ic regulated power supplies 0-5vdc150ma & 0- 20vdc150ma.

Three mr65 squaredigital meters for voltage & current measurement.

One mosfet irf540 mounted behind the panel.

Instrument fitted in highly shockproof plastics box.

Instrument dimensions: 12"x8"

accessories:

sufficient no of 4mm patch cords for interconnection

instruction manual for student workbook.

Power cords 230v, 50hz