

# Mosfet Characteristics Apparatus



**Product Categories:** [Characteristics And Application](#), [Electronics](#)

**Product Page:** <https://www.labappara.com/product/mosfet-characteristics-apparatus/>

## Product Description

Mosfet Characteristics Apparatus

**Features :**

Instrument comprises of Two DC Regulated Power Supplies 0-5 V DC/150mA & 0-25V DC/150 mA, three round meters for voltage & current measurement, one MOSFET IRF 840 mounted behind the panel, connections of Supplies, Meters & MOSFET brought out at 4mm Sockets.