


Multi Control Trainer



Product Categories: [Engineering Equipment](#), [Heat Transfer Lab](#)

Product Page: <https://www.labappara.com/product/multi-control-trainer/>

Product Description

Multi Control Trainer

Multi Process Control Trainer is the most flexible and modular system in its range. It is meant for studying various control loops used in industrial processes. Fabricated using high quality material, it can be used for elaborating the study of Controlling Different Process parameters like flow, level, pressure and temperature.

Features:

Self explanatory systems

Electrical safety consideration

Inbuilt safety measures

Innovative design

Feedback pressure control

Feedback temperature control

System frame with caster wheel arrangements

M.S powder coated cubical plant with standard instrument mountings

Specifications:

Electric Supply : 1 ϕ 230 V AC, 50 Hz

System Dimensions : 6 Ft. (L) X 2.5 Ft. (W) X 5.5 Ft. (H)**Applications :**

Study of PLC

Study of PLC based control of Multiprocess systems

Study of data logging system

Study of cascade, feedback, ratio and on-off, control schemes

Calibration of various sensors and transmitters

Study of SCADA based applications