

Multi Purpose Geo Sticks Type 2



Product Categories: [Geometry Items](#), [Math Lab](#)

Product Tags: [math lab instruments](#)

Product Page: <https://www.labappara.com/product/multi-purpose-geo-sticks-type-2/>

Product Description

Multi Purpose Geo Sticks Type 2 - Mathematics Laboratory Equipment

Moulded Plastic Strip having 2 holes of dia 5 mm at a distance of 50 mm and 120 mm from one end. 3 slots 125 x 5 mm with both ends rounded to semicircle. Location of slots(1st slot : 0-25 mm, 2nd slot : 70-95 mm, 3rd slot: 140-235 mm).