

## Newton Color Disc



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/newton-color-disc/>

### Product Description

#### **Newton Color Disc: Physics Lab Equipment**

Newton Color Disc for demonstrating that white light may be composed from all the spectral colors by rotation of multi color disc.

Comprising of a disc 200 mm, mounted on metal stand and driven by belt with the help of driving wheel fitted with a handle, can be clamped on table.

#### **Specifications:**

Disc 200 mm

Complete assembly mounted on a stable base.