


Newton's Ring Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/newtons-ring-microscope/>

Product Description

Newton's Ring Microscope - Scientific Equipment

Compact instrument with full arrangements for conducting Newton's Ring experiment.

A standard microscope unit with 30 x magnification and rotatable cross line eyepiece.

Focusing of microscope is through Rack and Pinion.

Longitudinal movement of 26 mm is done by rotating the drum provided which reads up to 0.001 cm.

Complete with optical glass and convex lens inside a metal case.

Packed in a sturdy wooden box along with all standard accessories.