

NICHOLSON HYDROMETER MEDIUM



Product Categories: [Heat And Thermodynamics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/nicholson-hydrometer-medium/>

Product Description

NICHOLSON HYDROMETER MEDIUM - Physics Equipment.

For Solid and liquid density experiments. It consists of a hollow metal cylinder, upper loading pan and bucket suspended of the bottom for stability and to hold the specimen, when immersed Made from brass sheet.

Available in 2 sizes, Medium & Large

Nicholson's hydrometer is a specialized device designed by our expert team for specific gravity, volume and weight related experiments.

It is a major device for the physical science experiments and requires precision as well as perfection.

This is the original hydrometer which is made up of brass core and finished extremely well by baked enamel finish.

This hydrometer is highly recommended for depth related experiments too and hence it finds a wide scale application.

Product data have been exported from - [Labappara scientific instrument](#)

Export date: Fri Jan 31 0:25:04 2025 / +0000 GMT