

Noddle Slide Assembly



Product Categories: [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/noddle-slide-assembly/>

Product Description

Noddle Slide Assembly - Physics Equipment

Noddle Slide Assembly Consists of a combination of Signal & double lens holders on a mount routing along the vertical axis. The combination lens holders are adjustable in height , and the single & twin lens holders can either be moved together by rack & pinion or relative to each other by sliding. These motions can be read by individual pointers on am scale.

(a) Plain, inalienable with 3" steel C.P Stand

(b) With fine machined cross slit, one surface white & the other optical black, with 3" steel C.P stem

c) An electric bulb is encased in a metal case with cord & pin It works 220 VAC.

Features:

Optimum strength

Lightweight

Perfect finish

Specifications:

Optical Bench: 1.5 meter