

## Nozzle Meter Test Rig



**Product Categories:** [Engineering Equipment](#), [Hydraulic Lab](#)

**Product Page:** <https://www.labappara.com/product/nozzle-meter-test-rig/>

### Product Description

#### Nozzle Meter Test Rig

Nozzle Meter Apparatus : The apparatus consists of a Nozzle meter. The flow can be regulated by the valve provided at the down stream Present set-up is self-contained water re-circulating unit, provided with a sump tank and a centrifugal pump.

Flow control valve and by-pass valve are fitted in water line to conduct the experiment on different flow rates. Flow rate of water is measured with the help of measuring tank and stop watch Experiments To measure discharge through Nozzle meter and to calibrate it Utilities Required Electric supply 0.5 kW, 220 V AC, Single Phase Water supply Tap water connection.