

OHM'S Law Apparatus



Product Categories: [Lab Meters](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/ohms-law-apparatus/>

Product Description

OHM'S Law Apparatus- For Physics Laboratory

The apparatus is used to study and verify the ohms law

The Instrument comprises of the following:

Two moving coil meters to measure voltage & current

One wire wound Potentiometer.

Two 4 mm socket terminals are provided for the connection of unknown resistance

Two 4 mm socket terminals are provided for DC supply

Circuit diagram Printed on the front panel

One ON/OFF Switch is provided on the front panel

To verify ohms law consists of one dial type

Terminal one each sensitive voltmeter 5 volt & millimeter 500 ma, mounted in

spray painted metal box with terminal.

With step down transformer giving volt and lead wire.

Available in volts - 1.5 v, 3, 4.5 volt