

Operational Amplifier As Instrumentation Amplifier



Product Categories: [Electronics](#), [Engineering Equipment](#), [Operational Amplifier's](#), [Operational Amplifier's](#)

Product Page:

<https://www.labappara.com/product/operational-amplifier-as-instrumentation-amplifier/>

Product Description

Operational Amplifier As Instrumentation Amplifier

Features :

Instrument comprises of One ± 15 VDC Regulated Power Supply, 2 DC Regulated short circuit and overload protected continuously variable Power supplies, 3 no. IC 741 placed inside & connections brought out at sockets. One Potentiometer and combination of Resistances placed on the front panel.