

## Operational Amplifier As Inverting & Non-Inverting Amplifier



**Product Categories:** [Electronics](#), [Engineering Equipment](#), [Operational Amplifier's](#), [Operational Amplifier's](#)

**Product Page:**

<https://www.labappara.com/product/operational-amplifier-as-inverting-non-inverting-amplifier/>

### Product Description

Operational Amplifier As Inverting & Non-Inverting Amplifier

#### **Features :**

Instrument comprises of DC Regulated Power Supplies of 0-1.5 VDC/150 mA,  $\pm 15V$  DC & two meters to measure input & output voltage, circuit diagram is printed.

Dimension 10"x7"x3".

Digital Multimeter 3.5 Digit, Digital Voltmeter 20 VDC.