

Operational Amplifier As Voltage and Current Feedback Amplifier



Product Categories: [Electronics](#), [Engineering Equipment](#), [Operational Amplifier's](#), [Operational Amplifier's](#)

Product Page:

<https://www.labappara.com/product/operational-amplifier-as-voltage-and-current-feedback-amplifier/>

Product Description

Operational Amplifier As Voltage and Current Feedback Amplifier

The instrument comprises of the following built in parts:-

- 1) Fixed DC Regulated Power Supply of +15V.
- 2) DC Regulated variable power supply of 0 -1.5V.
- 3) IC 741 is placed inside the cabinet and connections are brought out at sockets.
- 4) Two dual mode range moving coil meters (1.5v/1.5mA and 15v/15mA) to read the output voltage and current are mounted on front pannel.
- 5) Different types of resistances and a potentiometer are also provided on front pannel.