

Optical Bench Double Rod(Half shaper)



Product Categories: [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/optical-bench-double-rod-half-shaper/>

Product Description

Optical Bench Double Rod(Half shaper) - Physics Equipment

Optical Bench Consists of two metal rods one of the evenly divided in mms supported on heavy cast iron feet. Four heavy metal slider and the accessories. Two of the sliders have Transverse motion. One meter long consists of two rods on heavy sliders, two sliders have liter motion with two lens holder & two needles.

With Following Accessories :-

- (a) **Half Shaper**:-. 1 Meter steel C.P Rod.
- (b) **Half Shaper** :- 1(1/2) Meter steel C.P Rod.
- (c) **Half Shaper** :- 1 Meter with Both Rod Brass Pipe Covered.
- (d) **Half Shaper** :- 1(1/2) Meter With Both Rod Brass Pipe Covered