


# Optical Bench Single Rod



**Product Categories:** [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/optical-bench-single-rod/>

## Product Description

Optical Bench Single Rod - Physics Optical Equipment

Comprises a meter long, square steel rod 18 mm thick graduated in mm and supported on two sturdy metal feet leveling screw. Four slides and four accessories ,two of the slider have fine transverse motion. One meter long, single square steel rod accurately graduated fitted on two sturdy metal feet with leveling screws. Four special designed heavy sliders having fine index marks. Two of the sliders are having fine transverse slow motion arrangement. Complete with four accessories (Comprises a meter long, square steel rod 18 mm thick graduated in mm and supported on two sturdy metal feet leveling screw. Four slides and four accessories, two of the slider have fine transverse motion.)