

Optical Bench With Wooden Bed And Metal Accessories



Product Categories: [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/optical-bench-wooden-bed-metal-accessories/>

Product Description

Optical Bench With Wooden Bed And Metal Accessories - Physics Optical Instruments

Optical bench contains one meter long seasoned wood scale in mms on heavy sports four metal base slides and four metal accessories. Two of the top slides are provided with transverse slow motion arrangement. The optical bench is commonly used in physics labs today, and consists of a long, rigid member with a linear scale applied to it. Holders for light sources, lenses and screens are placed on the apparatus so that image formation can be observed.

Features:

Stable base

Robustness

Rust resistance