

Optical Lever (Single Aluminum Body)



Product Categories: [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page:

<https://www.labappara.com/product/optical-lever-single-aluminum-body/>

Product Description

Optical Lever (Single Aluminum Body) - Physics Optical Equipment

Optical lever is used for measuring the thickness of thin articles by optical method.

Optical work 40 mm die plain mirrors mounted on a metallic arc. The arc is mounted on aluminum metallic base with two fixed legs and one leveling screw. This instrument is used to determine thermal conductivity.