

Organ Bath



Product Categories: [Physiology Apparatus](#)

Product Tags: [Laboratory Equipment](#)

Product Page: <https://www.labappara.com/product/organ-bath/>

Product Description

Organ Bath - Laboratory Equipment

Organ bath are made of transparent acrylic sheet.

And then complete with uprights, organ & oxygen tubes, capillary tubing, and front lever. Chemically welded to make experiment.

Temperature by immersion heater controlled by a thermostat

Organ bath of chamber

This product is widely used in various laboratory and hospitals due to its reliable performance and longer service life.

Product Features:

Leak proof dovetailed joints

Sturdy metal frame

Bath size of 8" x 4" x 5"