

## Organ Bath



**Product Categories:** [Physiology Apparatus](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/organ-bath/>

### Product Description

Organ Bath - Laboratory Equipment

Organ bath are made of transparent acrylic sheet.

And then complete with uprights, organ & oxygen tubes, capillary tubing, and front lever. Chemically welded to make experiment.

Temperature by immersion heater controlled by a thermostat

Organ bath of chamber

This product is widely used in various laboratory and hospitals due to its reliable performance and longer service life.

Product Features:

Leak proof dovetailed joints

Sturdy metal frame

Bath size of 8" x 4" x 5"