

## P.H Meter



**Product Categories:** [Pathology Lab](#), [Ph Meters](#)

**Product Tags:** [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/p-h-meter/>

### Product Description

P.H Meter - Laboratory Equipment

PH meter probes measure the voltage or potential difference of the solution in which it is dipped.

By measuring potential difference, hydrogen ion concentration can be calculated using the laboratory equation.

Which gives the relationship between Hydrogen ion concentration and Voltage or Potential.

#### **Feature -**

PH meter table model with electrode.

PH meter table model without electrode.

Robust design

Optimum quality

Easy to install

Flawless performance

Easy usage

Long service life

Specifications:

Display 3½ LED

Range 0-14 pH 0±1999 mV

Resolution 0.01 pH 1 mV

Temperature Compensation Manual Automatic

Power Requirement 220 V AC

With Dust Cover & Instruction Manual.