

## Parallel Operation Single Phase Transformer



**Product Categories:** [Electronics](#), [Engineering Equipment](#), [Experimental Setup For Electrical Lab](#), [Experimental Setup For Electrical Lab](#)

**Product Page:**

<https://www.labappara.com/product/parallel-operation-single-phase-transformer/>

### Product Description

Parallel Operation Single Phase Transformer

#### **Features :**

Exclusive and attractive designed panel

Stand alone operation

Designed by considering all the safety precautions

Diagrammatic representation for the ease of connections

On board high quality meters

Provided with an extensive e-manual

Technical Specification :

Mains supply : 230 V AC  $\pm 10\%$ , 50Hz Transformers (2Nos.) Rating : 1k VA Primary

Voltage : 0 - 230 V Secondary Voltage : 0 - 200 - 230 V Meters Used Voltmeter (MI) :

500 V (2 Nos.) Ammeter (MI) : 10 A (2 Nos.) MCB (Single Phase) : 10 A