

Paramecium Zoology Model



Product Categories: [Educational equipments](#), [Educational Models](#), [Zoology Models & Botany Model](#)

Product Tags: [Scientific Educational Model](#)

Product Page: <https://www.labappara.com/product/paramecium-zoology-model/>

Product Description

Paramecium Zoology Model- Scientific Educational Model Paramecium Zoology Model-Paramecium zoology Model

One side of the model clearly shows the oral groove and their continuation into the gullet and the other side shows the details of the internal structure, mounted on base, with key card. Overall dimensions 410 x 250 x 95 mm. approx.