

## Pencil Jockey



**Product Categories:** [General Physics](#)

**Product Tags:** [electronics](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/pencil-jockey/>

### Product Description

Pencil Jockey - Physics Laboratory Equipment

We have excellence in manufacturing, exporting & supplying an exclusive range of Jockey pencil type & general items. Pencil Jockey Is Used in physics Experiments such as Wheatstone bridge and potentiometer work.

Pencil Jockey comes with insulated handle, plated brass contact, and a plated brass terminal with 4 mm. socket.

Available in pencil jockey with brass parts.