

Petri Dishes



Product Categories: [Dishes Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#)

Product Page: <https://www.labappara.com/product/petri-dishes/>

Product Description

Petri Dishes - Laboratory Glassware Equipment

Available in O.D. x Ht mm. (approx.) -

these can be repeatedly autoclaved. The ridges on the inner side of the cap prevent the base from sticking to the cap.

For that, we are stated as the highly called for manufacturer of Petri Dish.

Acknowledged For :-

Molded in Polypropylene Clear & unbreakable Can be repeatedly autoclaved Ridges on the inner side Prevent the base from sticking to the cap

Features:

Precise size

Optimum strength

Seamless finish

High accuracy

Reliable performance

Corrosion resistant finish

Longer life

(a) 50 x 17

(b) 62 x 17

(c) 80 x 17

(d) 150 x 20

(e) 200 x 20