

## pH Meter (Pocket)



**Product Categories:** [Laboratory Equipments](#), [Ph Meters](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/ph-meter-pocket/>

### Product Description

#### pH Meter (Pocket) - Scientific Equipment

The pH of the water in your tank mix can sometimes affect the efficacy of pesticides. Insecticides, in particular, have a tendency to break down (hydrolyse) rapidly in alkaline water. Water pH can vary, depending on the source, from 5.0 to 9.5. Neutral water has a pH of 7.0, while alkaline water is higher than 7.0. If your water pH is much higher than 8.0, you may want to consider using an acidifying agent such as vinegar to lower the pH in the tank. Many of the pH-sensitive pesticides.

PH-80

Fast and Accurate

Simple to Use

Water Resistant

Features:-

Range 0-14 pH

Accuracy  $\hat{A} \pm 0.2$  pH

Resolution 0.1

Automatic temperature compensation

Calibration via offset trimmer pH 7