

## Photo Cell ( Valve Type) Characteristic Apparatus



**Product Categories:** [Educational equipments](#), [Electronics](#)

**Product Tags:** [electronics](#)

**Product Page:**

<https://www.labappara.com/product/photo-cell-valve-type-characteristic-apparatus/>

### Product Description

Photo Cell ( Valve Type) Characteristic Apparatus - Electronics Laboratory Instrument

Selenium Photo cell mounted in metal box. Supplied with light source with 60W lamp. Photo cell Kit with DC Micro-ammeter of 1000mA & DC Milli-voltmeter of 500mV mounted behind the panel, connections of Supplies, Meters are brought out on front panel.

Photo cell mounted in box having window for injecting light with light source.

Photo cell with Schematic Diagram.