

# Photo Diode Characteristics Apparatus



**Product Categories:** [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:** <https://www.labappara.com/product/photo-diode-characteristics-apparatus/>

## Product Description

Photo Diode Characteristics Apparatus

Features :

Instrument comprises of DC Regulated Power Supply 0-3 VDC/150 mA, two round meters for voltage & current measurements, Photo Transistor mounted in a metal box connections of Supplies & Photo Transistor brought out at 4mm Sockets. Lamp holder with 100 W& 50 cm wooden plank are given with this apparatus.