

Physical Balance



Product Categories: [Balance & Scale](#), [Chemistry Equipments](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/physical-balance/>

Product Description

Physical Balance - Laboratory Equipment

Which is used to find the mass or weight of the given object is physical balance.

This physical balance device can correct the mass of an object into milligram.

The amplitude of oscillations will decrease when the beam is raised and it becomes to rest at a point which is called the resting point.

Feature -

Used to find the mass or weight of the given object is physical balance.

Knife edges on polished wooden case

Double door case