

Pipette Pumps



Product Categories: [Laboratory Plastic Ware](#), [Pippetes Stand](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/pipette-pumps/>

Product Description

Pipette Pumps - Scientific Laboratory Plastic ware Equipment

Pipette Pump: presents a fast release pi-petting device for precise pi-petting & zip-quick emptying. The Knurled thumb wheel can be rotated to draw up or dispense liquids while the releaser may be pushed in for zip-quick emptying. The moulded chuck has a silicon rubber collet inside to hold various pipette of standard sizes. It resists acids & alkalies and can be easily disassembled for through cleaning.

Particulars	Packing
2 ml	12 Pc
10 ml	12 Pc
25 ml	12 Pc