


Pipettes Transfer Volumetric (One Mark)



Product Categories: [Laboratory Glassware](#), [Pipettes](#) , [Volumetric Flasks And Measuring Cylinder](#)

Product Tags: [glassware](#)

Product Page:

<https://www.labappara.com/product/pipettes-transfor-volumetric-one-mark/>

Product Description

Pipettes Transfer Volumetric (One Mark) - Laboratory Glassware Equipment

A volumetric pipette or bulb pipette allows extremely accurate measurement (to four significant figures) a volume of solution.

These pipettes have a large bulb with a long narrow portion above with a single graduation mark as it is calibrated for a single volume (like a volumetric flask).

Features:

Soft-touch tip ejection

Completely autoclavable

Easy to use

Utilized for conducting different researches

Durable

Capacity - 5 ml, 10 ml, 20 ml, 25 ml, 50 ml