

Plank's Constant Apparatus



Product Categories: [Educational equipments](#), [Electronics](#)

Product Tags: [electronics](#)

Product Page: <https://www.labappara.com/product/planks-constant-apparatus/>

Product Description

Planck's Constant Apparatus - Electronics Laboratory Equipment

Planck's constant apparatus is specially designed for determination of Planck's constant and work Function of Materials by Photoelectric Effect. Moreover it is useful to verify Inverse square law of radiation using a photoelectric cell.

Planck's Constant Apparatus with build in power supply .

And constant workable on 230 volts AC and photo electric cell mounted with separate terminals.

With in provided a light source

And two meters wooden optical bench.

Features:

Operational fluency

Compact design

Highly durable