

Plastic Storage Box



Product Categories: [Laboratory Plastic Ware](#), [Slide Storage](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/plastic-storage-box/>

Product Description

Plastic Storage Box - Scientific Laboratory Plastic ware Equipment

These Storage Boxes, made of polypropylene, prove to be an excellent option for storage of different powdered reagents & other material in the laboratory. Square shape of these boxes helps in saving shelf space. Available in three different sizes these boxes have polypropylene caps, are autoclavable, provide good chemical resistance & heavy duty service.