

Platinum Loop Holder



Product Categories: [Chemistry Equipments](#), [General Chemistry Labware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/platinum-loop-holder/>

Product Description

Platinum Loop Holder - Chemistry Equipment

Features:

School type in test tube.

With Bakelite handles.

Platinum Loop Holder with Brass rod and chuck with Bakelite Handle.

9 inches long with ebonite handle

Adjustable nickellated head to take wires not coarser than 20 swg.