

## Pointer Eyepiece



**Product Categories:** [Eye Piece](#), [Microscope](#), [Microscope Assecories](#), [Optical Instruments](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/pointer-eye-piece-2/>

### Product Description

#### Pointer Eyepiece - Scientific Equipment

Pointer eyepieces are easily fit into the microscope. Spare pointer eyepiece 10 x or 15 x or 5 x accurately centered & fitted with fine pointer. Supplied in velvet case. The eyepiece pointer is a very useful microscope accessories.

fashioned easily for use on virtually any microscope position to use a microscope has run into the problem.

#### **Features:**

Rigid design

Perfect strength

Low maintenance