

Porous Plate Round



Product Categories: [Educational equipments](#), [Laboratory Porcelain Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/porous-plate-round/>

Product Description

Porous Plate Round - Scientific Laboratory Equipment

Desiccator 1 plate:-

Perforated porcelain desiccator plate with a large hole in the center for easy placement or removal from a glass desiccator. Plate is glazed on one side.

Due to its inherent properties, Porcelain Lab ware is widely used in the chemistry laboratory. It has excellent resistance to all reagents except hydroelectric acid, and is highly resistant to thermal shock while having a very low coefficient of expansion.

Size:- 100 mm, 125 mm, 150 mm, 200 mm, 250 mm, 300 mm

To prolong the life of porcelain products:

Follow gradual heating/cooling rates

Use an oven or hot plate and an intermediate step when quicker heating/cooling rates are required

Avoid contact of heated ware with a cold surface

Carefully inspect lab ware prior to each use and do not use any product that appears defective.