

## Post Office Box



**Product Categories:** [General Physics](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/post-office-box-2/>

### Product Description

#### Post Office Box

Post Office Box - Physics Laboratory Equipment

Post office box is education quality & Standard quality.

Brass block, plug type, having 3 pairs of ratio coils, 10, 100, 1000 ohms and 16 series coil constant wire.

All plugs are interchangeable with accuracy 0.1%

Constant coil standard teak box  $\frac{1}{2}$ " 3 ratio.

And with managing coils  $\frac{1}{2}$ " blocks.

Most superior brass block  $\frac{9}{16}$ " manganian coils.

Available in coils -

With 3 ratio coils.

With 4 ratio coils