

## Post Office Box Dial Type



**Product Categories:** [General Physics](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/post-office-box-dial-type/>

### Product Description

Post Office Box Dial Type - Physics Laboratory Instrument

Post office box is constant coils.

6 dial each consisting of 10 steps.

Provided with four ratios of 1, 10, 100, 1000 ohms and 2 ratios dials.

And then with Manganian coils 5 dials.

Plug type, in polished wooden box with hinged lid and panel to show circuit connections.

Coils of constantan wire double silk covered are wound non- inductively. Three pairs of ratio arms are each 10, 100, 1000 ohms and resistance arms contains 16 coils from 1 to 5000 ohms with infinity plug.

Total value 11, 110 ohms.

Two spring keys having contact for battery and galvanometer fitted.

Accuracy of coils 100 to 500 ohms + 0.05% and of the rest coils + 0.1%.

Available in coils -

4 dials , Constant coils

4 dials , Manganian coils

**Features:**

Precise design

Simple usage

Resistant against rus