

Post Office Box



Product Categories: [Electricity, Measurement & Circuits](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/post-office-box/>

Product Description

Post Office Box - Physics Equipment

Post office box in export quality.

And then with brass block , plug key, having 3 pairs of ratio coils - 10, 100, 1000 and 16 series constant wire.

All plug are interchangeable with accuracy 0.1%

Constant coil standard teak box $\frac{1}{2}$ " 3 ratio.

And with managing coils $\frac{1}{2}$ " blocks.

Most superior brass block $\frac{9}{16}$ " manganian coils.

With 3 ratio coils.

with 4 ratio coils.