

Potentiometers



Product Categories: [Electricity](#), [Measurement & Circuits](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/potentiometers/>

Product Description

Potentiometers For Physics Lab.

Potentiometers on nicely polished plywood base.

Potentiometers with full sun-mica top fitted with soft wood.

And with graduated scale 100 cm x 1 mm with brass CP connecting strips and brass terminals.

Potentiometers with Pencil jockey

And then with spring action Aluminum Jockey.

Available in wire - 4 wire , 10 wire